

DMI Medical Garments


Medical Compression Garments

Measuring Guidelines - Male

Measurements for Men's Compression Garments:

In order to help us get the best possible fit for the garment, please take the following measurements. Please take the measurements in **inches** and measure as indicated below.

Chest:	With arms relaxed at sides, measure around the chest going across the nipples and under the arm pits.
Waist:	Measure around the body going across the belly button.
Hips:	Measure around the fullest part of the lower buttocks so when the tape comes around to the front it is just above the pubic area.
Note:	To ensure a proper fit, please have someone else take your measurements and avoid wearing bulky clothes. Your feet should be 2-3" apart and your arms should be relaxed at your sides. We recommend the use of a flexible cloth measuring tape. Hold the tape level and snug but do not allow it to "dig" into your skin.



Custom Garments: If you do not find the size you are looking for, DMI can produce a custom garment for you. Please call us at 888-364-8721 and speak with a Customer Service Representative. You can also email us at: Info@DmiMedicalGarments.com